Pulse Oximeter

1byone Products Inc.
1230 E Belmont Street, Ontario, CA 91761
Customer Service: +1 909-391-3888
www.1byone.com

Instruction Manual/Bedienungsanleitung
Manuel d'instructions/Manual de instrucciones
Manuale di istruzioni

EN INTRODUCTION

DE Thank you for purchasing 1byone Pulse
FR Oximeter. This instruction manual
ES contains important information about
IT safety, usage, and disposal. Use the
product as described and keep this
manual for future use.

Precautions

1. Do not attempt to fix a broken Oximeter yourself. Only professionals with maintenance qualifications are allowed to perform interior maintenance as necessary.
2. Check the to-be-measured finger for skin integrity and proper blood circulation. Adjust the position of the sensor so the sensor is directly below the pad of the finger.

3. This product is not applicable to newborn babies.

4. Seek medical attention if the result of the measurement falls outside a normal healthy range and is not due to instrument malfunction.

5. Do not directly expose your eyes to the light-emitting components of the Oximeter as this can cause eye damage.

6. Please consult relevant medical literature carefully for details about clinical limitations and contraindications.

EN The following factors may influence the accuracy of the results:

- DE**
- FR**
- ES**
- IT**
1. Using the Oximeter near high-frequency devices, such as high-frequency electric knives and CT apparatuses.
 2. The Oximeter's probe is placed on the same arm as a blood pressure cuff, arterial duct, or intravenous injection.
 3. Users who suffer from hypotension, severe vascular atrophy, severe anemia, or hypoxemia.
 4. The user is suffering from a sudden cardiac arrest or in a state of shock.
 5. The user has applied nail polish or is using artificial nails.

Warnings

Do not use the Oximeter in an environment with any inflammable gases, inflammable anesthetic, or other inflammable substances.

Do not recharge used batteries as this can cause leakage, fire, or even explosion.

Do not use the Oximeter near an MRI or CT machine.

Do not operate the Oximeter when it is damp with water vapor condensation. Avoid moving the Oximeter from an excessively-cold environment to a high-temperature moist environment.

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Symbol	Description
	BF-type application part
	Caution: Please see this manual
%SpO2	Oxygen saturation
bpmPR	Pulse rate
	Temperature limitation
	When disposing of this product, send to an official recycling facility.

Overview

Oxygen saturation (SpO₂) is the percentage oxyhemoglobin (HbO₂) of all hemoglobin (Hb) and is an important physiological parameter in respiration, circulation, and the human body in general. The oxygen saturation of arterial blood is normally 98%. In general, the normal values of oxygen saturation should not be lower than 94%. If the measured value of oxygen saturation is lower than 94%, an insufficient amount of oxygen is being delivered to the body.

Pulse rate is a measure of heart beats per minute. The pulse rate of an adult is normally between 60 to 90 beats per minute.

The Perfusion Index (PI) is the ratio of pulsatile blood flow vs static blood in peripheral tissue. The PI of a normal adult is 3% or higher. Examination can

EN
DE
FR
ES
IT

EN still be performed even in a low or weak
DE PI condition.

FR

ES **Working Principles, Expected** IT **Usage, and Applicable Scope**

Based on full digital technology, the Pulse Oximeter non-invasively measures the actual content (oxygen saturation) of oxyhemoglobin (HbO₂) in arterial blood using the optical transmittance method.

The Pulse Oximeter measures the blood oxygen saturation and pulse rate of a human body via finger artery. It is applicable to a wide range of fields, such as home use, hospitals (including operation rooms in the departments of internal medicine and surgery, the department of anesthesiology, the department of paediatrics, and intensive care rooms), oxygen bars, social medical care institutions, and

centers for sports and health. Use this instrument for measurement before or after sports but it is advised not to use it during sports activities. Do not use this Oximeter for the long term measurement of patients.

EN

DE


FR

ES

IT

Schematic Diagram of Display

The following diagram shows the information display on the LED screen of the Oximeter while in a normal detection state:

- ❶ Oxygen saturation symbol
- ❷ Measured value of oxygen saturation
- ❸ Pulse amplitude bar graph
- ❹ SpO₂% plethysmograph signal
- ❺ Measured value of pulse rate
- ❻ Perfusion Index (PI) %
- ❼ Pulse rate symbol
- ❽ Battery power indicator

Power-On/Function Button

Press the Power-On/Function button to turn on the Oximeter. Press and hold this button again for one second to show the settings menus. "*" shows the currently selected setting. When "*" is at the top, press and hold to switch between the two setting screens (Alarm and Sound settings). Quick press to switch between the different setting options and press and hold on the selected option ("*") to change it.

Alarm and Beep Sound Setting

Hold the Power-On/Function button to access the two settings menus. These menus are shown below.

Alarm setup	*
Alarm	on
Beep	off
Demo	off
Restore	OK
Brightness	4
Exit	

Alarm Settings

Sounds Setup	*
SpO ₂ Alm Hi	100
SpO ₂ Alm Lo	94
PR Alm Hi	130
PR Alm Lo	50
+/-	+
Exit	

Sound Settings

Short press to move "*" to a setting option and then long press to change that option. When Alarm (Alm) is set to "on" and the measured values of blood oxygen saturation and pulse rate go beyond the set upper or lower limit, the Oximeter will play an alarm sound. When Alm is set to "off" and the measured values go beyond the set limit, the Oximeter will not play an alarm sound. When Beep is set to "on," a beep sound

EN will be played for every pulse beat during
DE pulse rate measurement. When Beep is
FR set to "off," no sound will be played
ES during pulse rate measurement.
IT

Note: If both Alarm and Beep are set to "on," only the beep sound of the pulse beats will be played. To make the alarm sound work properly, Beep must be set to "off."

Restore Factory Settings

In the Alarm Settings page, short press the Power-On/Function button to move "*" to Restore, then long press the button to restore factory settings. This setting will then read "OK."

Brightness Settings

In the Alarm Settings page, short press the Power-On/Function button to move "*" to Brightness, then long press the button to change between five

brightness settings, with 1 being the darkest and 5 being the brightest setting.

Alarm Range Settings


In the Sound Settings page, you can set the upper and lower limits of the SpO₂ Alarm and PR (Pulse Rate) Alarm. First, short press the Power-On/Function button to move "*" down to the +/- option so you can switch between increasing (+) or decreasing (-) the limits. Next, move "*" to each limit and increase or decrease the limit by long pressing the Power-On/Function button. Move "*" to the Exit option and long press to return to the monitoring interface.

Operation Guide

Switching on the Pulse Oximeter

Place a finger on the rubber pad inside the Oximeter. The pad of your finger should rest over the sensor and your

EN the fingernail (and the Oximeter's
 DE screen) should face upwards. Release
 FR the Oximeter so it closes on the finger.
 ES Finally, press the Power-On/Function
 IT button.

⚠ Measurements may be inaccurate if the pad of your finger is not placed entirely over the sensor.

Do not move your finger or body while measurement is ongoing.

The Oximeter will automatically turn off 10 seconds after it has been taken off your finger.

Display

While viewing the monitoring screen, press the Power-On/Function button to cycle through the 6 available screen displays. Display options are available for left and right hand use, as well as vertical viewing. These are shown below.

Replace the batteries when battery capacity is insufficient and the battery symbol () flashes on the screen.

Cleaning

Power off the Oximeter and remove the batteries before cleaning. Ensure the Oximeter is free of dust and dirt. Clean its surface, including LED screen, using 75% medical alcohol and a soft dry cloth.

Caution: Avoid liquid entering the Oximeter's body during cleaning.

Caution: Never immerse the Oximeter in liquid.

Disinfection

Wipe the Oximeter's rubber pads with a dry soft cloth dipped in 75% medical alcohol. The finger to be measured should also be disinfected using medical alcohol before and after measurement.

Do not disinfect the Oximeter by high-temperature, high-pressure, or gas disinfection.

Maintenance

- If the Oximeter will not be used for an extended period of time, remove the batteries from the battery compartment.
- Store the Oximeter in temperatures between 14° to 122°F (-10° to 50°C) and humidity no higher than 93%.
- Periodically check the Oximeter for damage.
- Avoid using the Oximeter in an environment with inflammable gases or using it in an environment where the temperature or humidity is excessively high or low.
- Check the accuracy of oxygen saturation and pulse rate readings by using an appropriate calibration apparatus.

Technical Specifications

1. Dimensions: 2.4x1.5x1.3in
(62.2x37x33.1mm)
Weight: 1.5oz (42.5g)
2. Peak wavelength range of the light emitted from the sensor: red light 660 nm \pm 3; infrared light 905 nm \pm 5.
3. Maximum optical power output of the sensor: 1.2 mW for infrared light (905 nm).
4. Normal working conditions:

Working temperature	41°F to 104°F (5°C to 40°C)
Relative humidity	15% to 80%, non-condensing
Atmospheric pressure	70 kPa to 106 kPa
Rated voltage	DC 3.0 V

5. Default values for alerts:

Parameter	Value
Oxygen saturation	Upper limit: 100 Lower limit: 94
Pulse rate	Upper limit: 130 Lower limit: 50
Alert condition	When the alert is turned on and the actual measured value goes beyond the preset alert parameter range, the Oximeter will play an alarm sound.


EN 6. Technical parameters:

Parameter		Value
Display range	Oxygen saturation	35% to 100%
	Pulse rate	25 bpm to 250 bpm
Resolution	Oxygen saturation	1%
	Pulse rate	1 bpm
Measurement precision	Oxygen saturation	±2% (70% to 100%) No requirement (≤ 69%)
	Pulse rate	±2 bpm
Alert range	Oxygen saturation	Upper limit: 50% to 100% Lower limit: 50% to 100%
	Pulse rate	Upper limit: 25 bpm to 250 bpm Lower limit: 25 bpm to 250 bpm
Alert error	Oxygen saturation	± 1% of the preset value
	Pulse rate	The greater of ±10% of the preset value and ±5 bpm
PI	Weak PI	Minimum 0.2%

EN Safety Standards

Type of anti-electric shock: internal power supply device
 Degree of anti-electric shock: BF-type application part
 Running mode: continuous working
 Waterproof grade: IP22

Storage and Transportation

EC REP Wellkang Ltd.
 Suite B, 29 Harley Street, London, W1G 9QR, UK

EN Disposal**DE Disposal of the appliance**

FR  Under no circumstances should you dispose of this product in normal domestic waste. This product is subject to the provisions of European Directive 2012/19/EU.

Disposal of batteries

Used batteries may not be disposed of in household waste.

All consumers are statutorily required to dispose of batteries at a collection point provided by their local municipality or retail store. The purpose of this obligation is to ensure batteries are disposed of in a non-polluting manner. Only dispose of batteries when they are fully discharged.

This product's packaging is made from environmentally friendly material and can be disposed of at your local recycling plant.

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