BODY ANALYSER SCALE

1byone Products Inc.
2313 E Philadelphia Street, Unit M, Ontario, CA 91761
www.1byone.com

Instruction Manual/Bedienungsanleitung
Manuel d'instructions/Manual de instrucciones/Manuale di istruzioni

1. Specification Parameters

Capacity: maximum 180kg/400lbs/28st 5lbs

Overload indication: Err

Low battery indication: Lo

Zero indication: 

High fat indication: ErrH

Low fat indication: ErrL

Number of user profiles: 10 users (P0-P9)

Age range: 10-80 years old

Height range: 100-240cm, 3'3"-7'10" ('ft,"in)

Impedance error: Err2

Measuring body fat range: 5%~50%

Measuring water range: 10%~80%


Measuring muscle range: 10%~80%

Measuring bone range: 0.5Kg~8.0Kg

Automatic identification range: 2Kg

Optional 'Athlete' mode: 


Displays BMI

2. Unit Conversion


2-1. When powering on the scale, the default measuring unit will be kg.

2-2. Users can choose kg, lbs, or st (stone). Press the SET button to wake up the scale, then press the UNIT button which is on the back of the scale to change the unit of measurement.


3. Set User Information


3-1. Press SET to enter settings mode.

3-2. Press UP or DOWN to change the user number and press SET to confirm.


3-3. Press UP or DOWN to change the Gender and Athlete modes. The symbol  means male and  means female. Press SET to confirm.


The symbol  means Athlete mode is selected. Athlete mode is for users who do physical at least 10 hours per week.

3-4. Press UP and DOWN to change age. Press SET to confirm.


3-5. Press UP and DOWN to change height. Press SET to confirm.


4. Weighing and measuring fat

4-1. After finishing confirming the settings, the scale will display "0.0".


4-2. Stand on the scale and wait for the weight measurement to settle on a final reading.

4-3. Next, four "0" will fade out one by one from left to right as the scale finalizes its body measurement calculations.


4-4. Once measurement is complete, the scale's display will cycle through the calculated BMI, body fat, body water, calories, bone mass, and muscle mass readings, finally displaying weight and user info again before automatically turning off.


5. Automatic identification

5-1. If a user's readings were correctly measured previously, the scale will remember the user's weight and automatically switch the scale to that user's profile the next time they use the scale.


5-2. Stand on the scale, wait for the weight measurement to settle on a final reading, and the scale will automatically switch to the correct user profile.


5-3. Next, four "0" will fade out one by one from left to right as the scale finalizes its body measurement calculations.

5-4. Once measurement is complete, the scale's display will cycle through the calculated BMI, body fat, body water, calories, bone mass, muscle mass, and weight.

5-5. If the weight of two or more users is very similar, the scale will display these user profiles (such as P1 and P2 in the sample graphic below). Press Up or Down to choose the correct user profile.

6. Evaluating body fat, body water, muscle mass, bone mass, and calorie readings for UW (Underweight) and OW (Overweight) users.

Chart 6-1

Body fat analysis (Unit: %)								
Age	Female				Male			
	UW	Healthy	OW	Obese	UW	Healthy	OW	Obese
10-39	<21	21-33	33-39	>39	<8	8-20	20-25	>25
40-59	<23	23-34	34-40	>40	<11	11-22	22-28	>28
60-80	<24	24-36	36-42	>42	<13	13-25	25-30	>30

Chart 6-2

Body water analysis (Unit: %)						
Age	Female			Male		
	UW	Healthy	OW	UW	Healthy	OW
10-15	<57	57-67	>67	<58	58-72	>72
16-30	<47	47-57	>57	<53	53-67	>67
31-60	<42	42-52	>52	<47	47-61	>61
61-80	<37	37-47	>42	<42	42-56	>56

Chart 6-3

BMI parameters (WHO)				
	Underweight	Healthy	Overweight	Obese
BMI	<18.5	18.5~25	25~30	>30

Chart 6-4

Body muscle mass analysis (Unit: %)		
Gender	Female	Male
Normal	>34%	>40%

Chart 6-5

Body bone mass analysis (Unit: Kg)						
Gender	Female			Male		
Weight	<45Kg	45-60Kg	>60Kg	<60Kg	60-75Kg	>75Kg
Bone	1.8Kg	2.2Kg	2.5Kg	2.5Kg	2.9Kg	3.2Kg

Chart 6-6

Calorie analysis (Unit: %)					
	Age	<17	18~29	30~50	>50
Male	<50Kg	1350	1200	1115	1075
	50~65Kg	1485	1320	1227	1183
	65~80Kg	1890	1680	1561	1505
	80~95Kg	2295	2040	1896	1828
	>95Kg	2700	2400	2230	2150
Female	<50Kg	1265	1180	1085	1035
	50~65Kg	1392	1298	1194	1139
	65~80Kg	1771	1652	1519	1449
	80~95Kg	2151	2006	1845	1760
	>95Kg	2530	2360	2170	2070

Disposal

Disposal of the appliance

Under no circumstances should you dispose of the appliance in normal domestic waste. This product is subject to the provisions of European Directive 2012/19/EU.

Dispose of the appliance via an approved disposal company or your municipal waste facility. Please observe currently applicable regulations. Please contact your waste disposal center if you need any further information.

Disposal of batteries

Used batteries may not be disposed of in household waste.

All consumers are statutorily required to dispose of batteries at a collection point provided by their local municipality or retail store. The purpose of this obligation is to ensure batteries are disposed of in a non-polluting manner. Only dispose of batteries when they are fully discharged.

The appliance's packaging is made from environmentally friendly material and can be disposed of at your local recycling plant.

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